

## Compact inspection system; profitability, safety, tracability



(Cabinet is optional)

The VizenDE is the compact version of the VIZEN, an automatic inspection system, which is essential within automated medication dispensing systems.

Are you still checking medication sachets manually? It is generally accepted that automation within medication distribution reduces the chance of errors. With the VizenDE, labour-intensive manual inspection is a thing of the past. The VizenDE ensures increased patient safety and improved efficiency in the production process.

The VizenDE checks the contents of the medication sachets based on photos the system takes using a powerful dome light. The medication in the sachets is assessed on quantity, size, colour and shape, after which the contents are compared to the order from the software APMC.

As the information on the inspected medication sachets can be called up at any time, each sachet is subsequently fully traceable.

## VizenDE


Automatic Medication Inspection



# VizenDE



**Quick and accurate medication inspection system**  
 The VizenDE automatically inspects and compares all medication in the medication sachets with information in the order. This comparison is based on quantity, size, colour and shape. All information is registered as an image and stored in the database to support safe and reliable administration of medication.



**Compact system**  
 As the VizenDE is compact, it's easy to integrate into your workspace. The VizenDE fits on a table for example. The VizenDE can also be delivered with matching cabinet.



**Management of the inspection results**  
 It is possible to assess rejected sachets based on the information and photo on the VizenDE. The information and the photos of the medication sachets are stored.



**Adjustable monitor**  
 The six-way, adjustable screen improves user-friendliness. The screen is large enough to view the images properly.



**Status and error reports**  
 The VizenDE is controlled through a 17" touch-screen monitor. Inspection reports and system alerts are displayed on the monitor and indicated on a signal lamp, so you are always aware of the status of the VizenDE.



**Registration**  
 The registration of new medications can be carried out on the VizenDE without use of an external system.

## Features & Benefits

- Quick inspection
- Generates detailed inspection reports
- Automatic investigation of suspected errors
- Control of an average of 40 medication sachets per minute
- Profitable
- Compact system
- Optimises the process of inspections
- Is integrated with APMC software, no extra computer is needed

## Specifications

VizenDE		
	<b>Model</b>	JV-VD10
	Dimensions (w x d x h)	608 x 486 x 540 mm
	Weight	around 63 kg - maximum of 75 kg depending on options
	Display	Touch screen 17"
	Electricity consumption	120W, max. 250W
	Power supply	100 ~ 240V -, 50/60Hz
	Recommended temperature	0 - 40 °C
	Humidity	10 ~ 80 %
	Optional	Matching cabinet

The specified values are target values. Right to deviations is retained. No rights can be derived from the technical specifications.